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Tunstall

Intended use

The connection sockets described in these installation instructions for installation in medical supply units is intended for the connection to the RAN room bus with speech in a nurse call system from the Flamenco system family. The connection sockets have a call button as well as a magnetic socket and an 8-pin socket for connecting operating devices.

Table 1. Products

70 0437 00	Connection socket combi, bedhead unit MAG
70 0437 50	Connection socket combi, bedhead unit, TVL MAG

Safety notes

NOTICE

The complete installation of the system is described in the technical manual.

CAUTION

The PCB is fitted with electrostatically sensitive components. Therefore, avoid direct contact.

In medical supply units, the provisions of DIN EN ISO 11197 apply to the laying of the nurse call system cables.

Product description

The **magnetic socket** is used to connect a device with a magnetic plug, e.g.

- ePAT[®]lite with magnetic plug (77 0370 10)
- ePAT[®]lite (77 0370 00) with magnetic adapter (77 0381 00)
- Patient handset with magnetic plug (76 0720 00)
- Pear push switch with magnetic plug (76 0710 0x, 76 0711 0x)


The **8-pin socket** is used to connect a device with an 8-pin mini DIN plug, e.g.

- Pear push switch (70 0710 0x, 70 0711 0x, 70 0712 00)

The **red call button** is used to raise calls. The call button lights up dimly to help you find the button in the dark. It lights up brightly as soon as a call is triggered.

The connection supports the following functions of the connected devices: transmission of calls, control of two light sources, TV transmission, Entertainment transmission, blinds control.

Figure 1. Control elements

- [1] Call button (red)
- [2] Magnetic socket
- [7] 8-pin socket

For a description of the remaining item numbers, see Figure 2.

Lighting control

WARNING

The use of inadmissible relays can endanger people's lives! Observe the following requirements for the light relays.


The switching outputs LT1 (reading light) and LT2 (room light) provide an output voltage of 24 V DC and a current of max. 40 mA, which is taken from the power supply of the room terminal. The potential refers to the 0 V connection of the room terminal. The switching pulse (24 V DC, 40 mA) supplied by the outputs lasts as long as the button on the control device (ePAT[®]lite, Patient handset, pear push switch) is pressed.

Table 2. Light relay requirements

Application	Switching light on/off (without dimming)	Dimming the light
Relay type	• Pulse relay (electronic) • Pulse relay (mechanical)	• Switching relay (electronic) • Switching relay (mechanical)
Nominal control voltage	24 V DC	
Control voltage range	18 – 26 V DC	
Maximum current consumption	40 mA	
Freewheeling diode	 When using mechanical relays, a freewheeling diode (e.g. 1N4002) directly on the relay is essential.	
Potential separation	Safe isolation with 2 x MOPP according to DIN EN 60601-1 must be provided.	
Ballast	-	For the dimming function, a corresponding dimmable ballast is required (e.g. OS-RAM DALI).

Setting jumpers

Figure 2. Jumpers and connections

- [3] Socket for flat cable for connection to the connection panel of the medical supply unit
- [4] Jumpers for TV sound channel
- [5] Jumper for blinds option
- [6] Jumper for light option
- [8] Socket for connection light control reading light, max. 40 mA. Required cable with plug not included: article no. 50 0308 00.
- [9] Jumpers for bed number

For a description of the remaining item numbers, see Figure 1.

Setting the bed number

Bed number	Jumpers inserted for bed number
Bed 1	P1
Bed 2	P2
Bed 3	P1, P2 (= factory setting)
Bed 4	P3
Bed 5	P1, P3
Bed 6	P2, P3

Setting the light option

Light option	Jumper for light option
The light buttons LT1 and LT2 of the operating device switch different light sources. The outputs LT1 (reading light) and LT2 (room light) are separate. Two light relays are required.	Jumper removed
Both light buttons (LT1 and LT2) of the operating device switch the same light source. Outputs LT1 and LT2 are connected. Only one light relay is required.	Jumper inserted

Setting blinds option

Blinds option	Jumpers for blinds option
No blinds control connected.	Jumper 1 inserted
Blinds control.	Jumper 2 inserted

NOTICE

The blinds control is not connected to the connection socket combi, bedhead unit MAG, but to a universal RAN interface (70 0848 00), which in turn is connected to the RAN room bus.

Setting the TV sound channel

Two jumpers must always be inserted to set the TV sound channel.

If there are several TVs in the room, set the sound channel for each bed that you want it to receive.

Bed number	Jumpers inserted for TV sound channel
Bed 1	TV1
Bed 2	TV2
Bed 3	TV3
Bed 4	TV4


If you want all beds to receive sound from the same TV set, set Sound Channel 1, i.e. Jumpers TV1.

Mounting

The installation of the connection socket in the medical supply unit depends on the particular design of the supply unit.

1. Open the housing of the medical supply unit.
2. Insert the connection socket.
3. Plug the flat cable into the socket for connection to the connection panel of the medical supply unit.
4. If there is a reading light, plug the cable (50 0308 00) plug into the socket for light control reading light connection.
5. Close the housing of the medical supply unit.


Figure 3. Light control connection reading light

Medical supply unit connection panel

The connection socket must be connected to the connection panel of the medical supply unit via a 40-pin flat cable. The connection points are assigned as follows:

Figure 4. Connection points of the connector for the flat cable

Notes on the connections

Wire bridges

- All connection points **+24V-Si** in the connection panel must have the same potential, i.e. either be connected internally or connected by external wire bridges.
- All connection points **0V** in the connection panel must have the same potential, i.e. either be connected internally or connected by external wire bridges.

RAN room bus

- The cable for speech transmission, i.e. the cable containing the wires **MIC**, **LS**, and **SCH-B**, may be stripped for a max. of 30 mm only, to ensure the protective screen function. Otherwise, a short circuit situation or crosstalk may occur.
- The cable shield and the shield wire of the cable containing the wires **RAN**, **+24V-Si** and **0V** are not connected and must be removed in the back box up to the cable sheath.
- The unused wires must not be interconnected.

Connections

Figure 5. Connections to the connection panel of the medical supply unit